

UNRAVELING THE CONFORMATIONAL DYNAMICS OF STAPHYLOKINASE

Naina Verma^{1,2*}, Lenka Hernychova^{1,2}, Jiri Damborsky²

1) RECAMO, Masaryk Memorial Cancer Institute, Zluty kopec 7, 656 53 Brno, Czech Republic

2) Loschmidt Laboratories, Department of Experimental Biology and RECETOX, Faculty of Science, Masaryk University, Kamenice 5, 625 00 Brno, Czech Republic
Email: naina.verma@mou.cz

BACKGROUND

The structural proteomic approach, employing hydrogen-deuterium exchange coupled with mass spectrometry (HDX-MS), supported by bioinformatics and computational modelling can unravel a plethora of information about protein dynamics. Using this state-of-the-art technique in understanding the in-solution conformational dynamics of a potential thrombolytic drug, Staphylokinase (SAK), will play a crucial role in providing important targets for protein engineering to develop a highly efficient, and alternative thrombolytic drug to Alteplase, which remains the sole FDA-approved drug for the treatment of ischemic stroke.¹



OBJECTIVES

- To decipher the in-solution dynamics of SAK in:
 - Apo form (wild-types and mutants)
 - Holo form (wild-types and mutants with interaction partners: plasmin and plasminogen)
- To understand the impact of mutations on SAK dynamics.

METHODOLOGY

In-solution conformational dynamics of SAK

- Hydrogen-deuterium exchange coupled with mass spectrometry**
The structure of SAK in its unbound and bound forms with plasmin and plasminogen will be determined.
- This will provide baseline data to further improve the design of mutants.



Protein engineering and biochemical characterization

- Computational and experimental approaches:** Afflib, Ribosome display, machine learning, RF diffusion
- Computational validation** by MD simulations
- Fibrin plate assay**: clot dissolving activity
- Pharmacokinetics**: to check half-life
- Immunogenicity testing**: *in-vitro* assays

EXPERIMENTAL DESIGN

Peptide mapping


- primary sequence of the protein
- sequence coverage of the sample
- peptide list for HDX data analysis

H/D exchange


- amide H⁺ replaced by D
- depends on the strength of H-bond
- reveals the conformational dynamics

Data evaluation

- HDExaminer software
- calculation of single amino acid deuterium uptake in peptides
- tools to visualize evaluated data



PRELIMINARY RESULTS



These 3D models of SAK show the variations in deuterium uptake (**red-maximum, blue-minimum, and black-uncovered peptides**) of SAK mutants compared to the wild-type. (A) SAK SY155 exhibits overall decreased solvent accessibility despite higher deuteration observed in α -helix region while (B) SAK THR144 is more flexible with a significant increase in deuteration in the loop and exhibits higher solvent accessibility.

OUTCOMES

- In-solution conformational dynamics of SAK.
- In-solution conformational dynamics of SAK with interacting partners.
- Physiologically compatible conformation of SAK for enhanced and efficient thrombolysis.

REFERENCES:

1. Nikitin, D. *et al.* Computer-aided engineering of staphylokinase toward enhanced affinity and selectivity for plasmin. *Comput. Struct. Biotechnol. J.* 20, 1366–1377 (2022).

ACKNOWLEDGEMENTS:

This study is supported by MH CZ - DRO (MMCI, 00209805) and the Czech Science Foundation, project no. 22-09853S. The author(s) thank the funding agencies. Views and opinions expressed are however those of author(s) and do not necessarily reflect those of funding agencies and the European Union.